





RLBP1 Antibody

Product Code	CSB-PA019743LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12271
Immunogen	Recombinant Human Retinaldehyde-binding protein 1 protein (10-91AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Retinaldehyde-binding protein 1 (Cellular retinaldehyde-binding protein), RLBP1, CRALBP
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RLBP1

Image



Immunohistochemistry of paraffin-embedded human eye tissue using CSB-PA019743LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA019743LA01HU at dilution of 1:100 and Cy3-conjugated Goat Anti-Rabbit IgG



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.