



RLBP1 Antibody, FITC conjugated

Product Code	CSB-PA019743LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12271
Immunogen	Recombinant Human Retinaldehyde-binding protein 1 protein (10-91AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Retinaldehyde-binding protein 1 (Cellular retinaldehyde-binding protein), RLBP1, CRALBP
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RLBP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.