



# RLBP1 Antibody, FITC conjugated

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA019743LC01HU   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                              |
| <b>Uniprot No.</b>         | P12271   |
| <b>Immunogen</b>           | Recombinant Human Retinaldehyde-binding protein 1 protein (10-91AA)                        |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Human  |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | FITC   |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4           |
| <b>Purification Method</b> | >95%, Protein G purified   |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Alias</b>               | Retinaldehyde-binding protein 1 (Cellular retinaldehyde-binding protein),<br>RLBP1, CRALBP |
| <b>Immunogen Species</b>   | Homo sapiens (Human)   |
| <b>Research Area</b>       | Neuroscience   |
| <b>Target Names</b>        | RLBP1  |