





RPS21 Antibody

Product Code	CSB-PA020399LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WVC2
Immunogen	Recombinant Human 40S ribosomal protein S21 protein (1-81AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	40S ribosomal protein S21 (Small ribosomal subunit protein eS21), RPS21
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RPS21

Image



Immunohistochemistry of paraffin-embedded human endometrial cancer using CSB-PA020399LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA020399LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.