



SLC25A1 Antibody, Biotin conjugated

Product Code	CSB-PA021475LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53007
Immunogen	Recombinant Human Tricarboxylate transport protein, mitochondrial protein (144-182AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tricarboxylate transport protein, mitochondrial (Citrate transport protein) (CTP) (Solute carrier family 25 member 1) (Tricarboxylate carrier protein), SLC25A1, SLC20A3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC25A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.