




TXNL4A Antibody

Product Code	CSB-PA025386LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P83876
Immunogen	Recombinant Human Thioredoxin-like protein 4A protein (1-142AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thioredoxin-like protein 4A (DIM1 protein homolog) (Spliceosomal U5 snRNP-specific 15 kDa protein) (Thioredoxin-like U5 snRNP protein U5-15kD), TXNL4A, DIM1 TXNL4
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TXNL4A

Image



Immunofluorescent analysis of PC-3 cells using CSB-PA025386LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)