



TXNL4A Antibody, HRP conjugated

Product Code	CSB-PA025386LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P83876
Immunogen	Recombinant Human Thioredoxin-like protein 4A protein (1-142AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thioredoxin-like protein 4A (DIM1 protein homolog) (Spliceosomal U5 snRNP-specific 15 kDa protein) (Thioredoxin-like U5 snRNP protein U5-15kD), TXNL4A, DIM1 TXNL4
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TXNL4A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.