



VNN1 Antibody, Biotin conjugated

Product Code	CSB-PA025883LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95497
Immunogen	Recombinant Human Pantetheinase protein (97-197AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pantetheinase (EC 3.5.1.92) (Pantetheine hydrolase) (Tiff66) (Vascular non-inflammatory molecule 1) (Vanin-1), VNN1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	VNN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.