





ZNF331 Antibody

Product Code	CSB-PA026680LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQX6
Immunogen	Recombinant Human Zinc finger protein 331 protein (179-442AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 331 (C2H2-like zinc finger protein rearranged in thyroid adenomas) (Zinc finger protein 361) (Zinc finger protein 463), ZNF331, RITA ZNF361 ZNF463
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZNF331

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA026680LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA026680LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.