



ZNF331 Antibody, HRP conjugated

Product Code	CSB-PA026680LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQX6
Immunogen	Recombinant Human Zinc finger protein 331 protein (179-442AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 331 (C2H2-like zinc finger protein rearranged in thyroid adenomas) (Zinc finger protein 361) (Zinc finger protein 463), ZNF331, RITA ZNF361 ZNF463
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZNF331
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.