




ZNF695 Antibody

Product Code	CSB-PA026957LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IW36
Immunogen	Recombinant Human Zinc finger protein 695 protein (1-88AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 695 (Zinc finger protein SBZF3), ZNF695
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZNF695

Image



Western Blot

Positive WB detected in: Rat liver tissue, Mouse liver tissue, Mouse kidney tissue


All lanes: ZNF695 antibody at 3.4µg/ml

Secondary


Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 61, 21, 16, 11 kDa


Observed band size: 61 kDa



Immunohistochemistry of paraffin-embedded human bladder cancer using CSB-PA026957LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human cervical cancer using CSB-PA026957LA01HU at dilution of 1:100



Immunofluorescent analysis of Hela cells using CSB-PA026957LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.