



I79_001071 Antibody, Biotin conjugated

Product Code	CSB-PA20599D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	G3GTT2
Immunogen	Recombinant Cricetulus griseus C-C motif chemokine protein (25-143AA)
Raised In	Rabbit
Species Reactivity	Cricetulus griseus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-C motif chemokine, I79_001071
Immunogen Species	Cricetulus griseus (Chinese hamster) (Cricetulus barabensis griseus)
Research Area	Others
Target Names	I79_001071
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.