



IX Antibody, Biotin conjugated

Product Code	CSB-PA301675LD01ECY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P69538
Immunogen	Recombinant Enterobacteria phage M13 Tail virion protein G9P protein (1-32aa)
Raised In	Rabbit
Species Reactivity	Enterobacteria phage M13
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tail virion protein G9P (Coat protein C, polypeptide II) (G9P), IX
Immunogen Species	Enterobacteria phage M13
Research Area	Others
Target Names	IX
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.