



IX Antibody, Biotin conjugated

Product Code	CSB-PA301675LD01ECY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P69538
Immunogen	Recombinant Enterobacteria phage M13 Tail virion protein G9P protein (1-32aa)
Raised In	Rabbit
Species Reactivity	Enterobacteria phage M13
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tail virion protein G9P (Coat protein C, polypeptide II) (G9P), IX
Immunogen Species	Enterobacteria phage M13
Research Area	Others
Target Names	IX