




HLA-DRB3 Antibody


Product Code	CSB-PA302920LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P79483
Immunogen	Recombinant Human HLA class II histocompatibility antigen, DR beta 3 chain protein (111-227AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:500-1:5000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	HLA class II histocompatibility antigen, DR beta 3 chain (MHC class II antigen DRB3), HLA-DRB3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	HLA-DRB3

Image



Western Blot

Positive WB detected in: Raji whole cell lysate,
 Mouse brain tissue
 All lanes: HLA-DRB3 antibody at 3.5µg/ml
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 30 kDa
 Observed band size: 30 kDa



Immunofluorescent analysis of A549 cells using CSB-PA302920LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.