



HLA-DRB3 Antibody, HRP conjugated

Product Code	CSB-PA302920LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P79483
Immunogen	Recombinant Human HLA class II histocompatibility antigen, DR beta 3 chain protein (111-227AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	HLA class II histocompatibility antigen, DR beta 3 chain (MHC class II antigen DRB3), HLA-DRB3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	HLA-DRB3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.