



afp Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA325933LC01APN |
| Abbreviation | Antifungal protein |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P17737 |
| Immunogen | Recombinant Aspergillus giganteus Antifungal protein (44-94AA) |
| Raised In | Rabbit |
| Species Reactivity | Aspergillus giganteus |
| Tested Applications | ELISA |
| Relevance | This protein inhibits the growth of a variety of fungal species. |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Antifungal protein, afp |
| Immunogen Species | Aspergillus giganteus |
| Research Area | Others |
| Target Names | afp |