




hupB Antibody

Product Code	CSB-PA363134HA01SWW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A1R9
Immunogen	Recombinant Salmonella typhi DNA-binding protein HU-beta protein (1-90AA)
Raised In	Rabbit
Species Reactivity	Salmonella typhi
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA-binding protein HU-beta (HU-1) (NS1), hupB
Immunogen Species	Salmonella typhi
Research Area	Others
Target Names	hupB

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: hupB antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 30 kDa

Observed band size: 30 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.