





TTC19 Antibody

Product Code	CSB-PA724947LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6DKK2
Immunogen	Recombinant Human Tetratricopeptide repeat protein 19, mitochondrial protein (108-380AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetratricopeptide repeat protein 19, mitochondrial (TPR repeat protein 19), TTC19
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TTC19

Image



Immunohistochemistry of paraffin-embedded human cervical cancer using CSB-PA724947LA01HU at dilution of 1:100



Immunofluorescent analysis of MCF-7 cells using CSB-PA724947LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.