



TTC19 Antibody, FITC conjugated

Product Code	CSB-PA724947LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6DKK2
Immunogen	Recombinant Human Tetratricopeptide repeat protein 19, mitochondrial protein (108-380AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetratricopeptide repeat protein 19, mitochondrial (TPR repeat protein 19), TTC19
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TTC19
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.