



ZFYVE26 Antibody, FITC conjugated

Product Code	CSB-PA734654LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q68DK2
Immunogen	Recombinant Human Zinc finger FYVE domain-containing protein 26 protein (2306-2473AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger FYVE domain-containing protein 26 (FYVE domain-containing centrosomal protein) (FYVE-CENT) (Spastizin), ZFYVE26, KIAA0321
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZFYVE26