




hist2h2l Antibody

Product Code	CSB-PA740767LA01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6PC60
Immunogen	Recombinant Zebrafish Histone H2B 3 protein (2-126AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Histone H2B 3, hist2h2l
Immunogen Species	Zebrafish
Research Area	Others
Target Names	hist2h2l

Image



Western Blot

Positive WB detected in: Zebrafish tissue 20µg, 10µg,

All lanes: hist2h2l antibody at 3.4µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 14 kDa

Observed band size: 14 kDa