



ZBTB33 Antibody, HRP conjugated

Product Code	CSB-PA769768LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86T24
Immunogen	Recombinant Human Transcriptional regulator Kaiso protein (113-228AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcriptional regulator Kaiso (Zinc finger and BTB domain-containing protein 33), ZBTB33, KAISO ZNF348
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZBTB33