



ZBTB33 Antibody, Biotin conjugated

Product Code	CSB-PA769768LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86T24
Immunogen	Recombinant Human Transcriptional regulator Kaiso protein (113-228AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcriptional regulator Kaiso (Zinc finger and BTB domain-containing protein 33), ZBTB33, KAISO ZNF348
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZBTB33
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.