



PXK Antibody, Biotin conjugated

Product Code	CSB-PA770363LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z7A4
Immunogen	Recombinant Human PX domain-containing protein kinase-like protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PX domain-containing protein kinase-like protein (Modulator of Na,K-ATPase) (MONaKA), PXX
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	PXK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.