



FAT3 Antibody, HRP conjugated

Product Code	CSB-PA819472LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDW7
Immunogen	Recombinant Human Protocadherin Fat 3 protein (1209-1470AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protocadherin Fat 3 (hFat3) (Cadherin family member 15) (FAT tumor suppressor homolog 3), FAT3, CDHF15 KIAA1989
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	FAT3