



FAT3 Antibody, HRP conjugated

Product Code	CSB-PA819472LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDW7
Immunogen	Recombinant Human Protocadherin Fat 3 protein (1209-1470AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protocadherin Fat 3 (hFat3) (Cadherin family member 15) (FAT tumor suppressor homolog 3), FAT3, CDHF15 KIAA1989
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	FAT3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.