



PKD1L1 Antibody, FITC conjugated

Product Code	CSB-PA823467LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDX9
Immunogen	Recombinant Human Polycystic kidney disease protein 1-like 1 protein (921-1164AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polycystic kidney disease protein 1-like 1 (PC1-like 1 protein) (Polycystin-1L1), PKD1L1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PKD1L1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.