



RNASEH2C Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA837441LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8TDP1 |
| Immunogen | Recombinant Human Ribonuclease H2 subunit C protein (1-164AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Ribonuclease H2 subunit C (RNase H2 subunit C) (Aicardi-Goutieres syndrome 3 protein) (AGS3) (RNase H1 small subunit) (Ribonuclease HI subunit C), RNASEH2C, AYP1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | RNASEH2C |