



RNASEH2C Antibody, FITC conjugated

Product Code	CSB-PA837441LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDP1
Immunogen	Recombinant Human Ribonuclease H2 subunit C protein (1-164AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribonuclease H2 subunit C (RNase H2 subunit C) (Aicardi-Goutieres syndrome 3 protein) (AGS3) (RNase H1 small subunit) (Ribonuclease HI subunit C), RNASEH2C, AYP1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RNASEH2C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.