




ISOC1 Antibody

Product Code	CSB-PA839302LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96CN7
Immunogen	Recombinant Human Isochorismatase domain-containing protein 1 protein (155-298AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Isochorismatase domain-containing protein 1, ISOC1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ISOC1

Image



Western Blot

Positive WB detected in: SH-SY5Y whole cell lysate


All lanes: ISOC1 antibody at 3.5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 33 kDa

Observed band size: 33 kDa



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA839302LA01HU at dilution of 1:100