



SLC7A6 Antibody, FITC conjugated

Product Code	CSB-PA849761LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92536
Immunogen	Recombinant Human Y+L amino acid transporter 2 protein (473-515AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Y+L amino acid transporter 2 (Cationic amino acid transporter, y+ system) (Solute carrier family 7 member 6) (y(+)-L-type amino acid transporter 2) (Y+LAT2) (y+LAT-2), SLC7A6, KIAA0245
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC7A6