



# SLC7A6 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA849761LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q92536
<b>Immunogen</b>	Recombinant Human Y+L amino acid transporter 2 protein (473-515AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Y+L amino acid transporter 2 (Cationic amino acid transporter, y+ system) (Solute carrier family 7 member 6) (y(+)-L-type amino acid transporter 2) (Y+LAT2) (y+LAT-2), SLC7A6, KIAA0245
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Neuroscience
<b>Target Names</b>	SLC7A6
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.