



EML4 Antibody, Biotin conjugated

Product Code	CSB-PA862066LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96HE3
Immunogen	Recombinant Human EML4 protein (1-62AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Echinoderm microtubule-associated protein-like 4 (EMAP-4) (Restrictedly overexpressed proliferation-associated protein) (Ropp 120), EML4, C2orf2 EMAPL4
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	EML4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.