



DNAJC4 Antibody, Biotin conjugated

Product Code	CSB-PA865087LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NNZ3
Immunogen	Recombinant Human DnaJ homolog subfamily C member 4 protein (1-155AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DnaJ homolog subfamily C member 4 (DnaJ-like protein HSPF2) (Multiple endocrine neoplasia type 1 candidate protein number 18), DNAJC4, HSPF2 MCG18
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DNAJC4