



DNAJC4 Antibody, Biotin conjugated

Product Code	CSB-PA865087LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NNZ3
Immunogen	Recombinant Human DnaJ homolog subfamily C member 4 protein (1-155AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DnaJ homolog subfamily C member 4 (DnaJ-like protein HSPF2) (Multiple endocrine neoplasia type 1 candidate protein number 18), DNAJC4, HSPF2 MCG18
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DNAJC4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.