




ABHD10 Antibody

Product Code	CSB-PA865140LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NUJ1
Immunogen	Recombinant Human Mycophenolic acid acyl-glucuronide esterase, mitochondrial protein (113-214AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mycophenolic acid acyl-glucuronide esterase, mitochondrial (EC 3.1.1.93) (Alpha/beta hydrolase domain-containing protein 10) (Abhydrolase domain-containing protein 10), ABHD10
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ABHD10

Image



Immunohistochemistry analysis of human small intestine tissue using CSB-PA865140LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.