



ABHD10 Antibody, HRP conjugated

Product Code	CSB-PA865140LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NUJ1
Immunogen	Recombinant Human Mycophenolic acid acyl-glucuronide esterase, mitochondrial protein (113-214AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mycophenolic acid acyl-glucuronide esterase, mitochondrial (EC 3.1.1.93) (Alpha/beta hydrolase domain-containing protein 10) (Abhydrolase domain-containing protein 10), ABHD10
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ABHD10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.