



MSRB1 Antibody, HRP conjugated

Product Code	CSB-PA882163LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZV6
Immunogen	Recombinant Human Methionine-R-sulfoxide reductase B1 protein (1-116AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Methionine-R-sulfoxide reductase B1 (MsrB1) (EC 1.8.4.12) (EC 1.8.4.14) (Selenoprotein X) (SelX), MSRB1, SEPX1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MSRB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.