



MSRB1 Antibody, Biotin conjugated

Product Code	CSB-PA882163LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZV6
Immunogen	Recombinant Human Methionine-R-sulfoxide reductase B1 protein (1-116AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Methionine-R-sulfoxide reductase B1 (MsrB1) (EC 1.8.4.12) (EC 1.8.4.14) (Selenoprotein X) (SelX), MSRB1, SEPX1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MSRB1