



GTPBP4 Antibody, HRP conjugated

Product Code	CSB-PA887173LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZE4
Immunogen	Recombinant Human Nucleolar GTP-binding protein 1 protein (347-634AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleolar GTP-binding protein 1 (Chronic renal failure gene protein) (GTP-binding protein NGB), GTPBP4, CRFG NOG1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GTPBP4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.