



WSB1 Antibody, Biotin conjugated

Product Code	CSB-PA897589LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6I7
Immunogen	Recombinant Human WD repeat and SOCS box-containing protein 1 protein (201-421AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	WD repeat and SOCS box-containing protein 1 (WSB-1) (SOCS box-containing WD protein SWiP-1), WSB1, SWIP1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	WSB1