





CTNND1 Antibody

Product Code	CSB-PA006172LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60716
Immunogen	Recombinant Human Catenin delta-1 protein (482-605AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Catenin delta-1 (Cadherin-associated Src substrate) (CAS) (p120 catenin) (p120(ctn)) (p120(cas)), CTNND1, KIAA0384
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CTNND1


Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA006172LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA006172LA01HU at dilution of 1:100



Immunofluorescent analysis of A549 cells using CSB-PA006172LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)