





CTSZ Antibody

Product Code	CSB-PA006206LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBR2
Immunogen	Recombinant Human Cathepsin Z protein (127-287AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cathepsin Z (EC 3.4.18.1) (Cathepsin P) (Cathepsin X), CTSZ
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	CTSZ

Image



Immunohistochemistry of paraffin-embedded human spleen tissue using CSB-PA006206LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA006206LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.