




ACTN4 Antibody

Product Code	CSB-PA00819A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43707
Immunogen	Recombinant Human Alpha-actinin-4 protein (6-277AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-actinin-4 (Non-muscle alpha-actinin 4), ACTN4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ACTN4

Image



Western Blot

Positive WB detected in: Mouse kidney tissue


All lanes: ACTN4 antibody at 3.4µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 105, 80, 60 kDa

Observed band size: 105, 80 kDa



Immunohistochemistry of paraffin-embedded human lung tissue using CSB-PA00819A0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.