





PTTG1IP Antibody

Product Code	CSB-PA00859A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53801
Immunogen	Recombinant Human Pituitary tumor-transforming gene 1 protein-interacting protein (118-180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pituitary tumor-transforming gene 1 protein-interacting protein (Pituitary tumor-transforming gene protein-binding factor) (PBF) (PTTG-binding factor), PTTG1IP, C21orf1 C21orf3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PTTG1IP

Image



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA00859A0Rb at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA00859A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.