



GCLM Antibody, Biotin conjugated

Product Code	CSB-PA009322LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48507
Immunogen	Recombinant Human Glutamate--cysteine ligase regulatory subunit protein (1-274AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutamate--cysteine ligase regulatory subunit (GCS light chain) (Gamma-ECS regulatory subunit) (Gamma-glutamylcysteine synthetase regulatory subunit) (Glutamate--cysteine ligase modifier subunit), GCLM, GLCLR
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GCLM
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.