



HAAO Antibody, Biotin conjugated

Product Code	CSB-PA010112LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46952
Immunogen	Recombinant Human 3-hydroxyanthranilate 3,4-dioxygenase protein (153-223AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	3-hydroxyanthranilate 3,4-dioxygenase (EC 1.13.11.6) (3-hydroxyanthranilate oxygenase) (3-HAO) (3-hydroxyanthranilic acid dioxygenase) (HAD), HAAO
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	HAAO
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.