



HSD17B7 Antibody, FITC conjugated

Product Code	CSB-PA010776LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56937
Immunogen	Recombinant Human 3-keto-steroid reductase protein (128-198AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	3-keto-steroid reductase (EC 1.1.1.270) (17-beta-hydroxysteroid dehydrogenase 7) (17-beta-HSD 7) (Estradiol 17-beta-dehydrogenase 7) (EC 1.1.1.62) (Short chain dehydrogenase/reductase family 37C member 1), HSD17B7, SDR37C1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HSD17B7