



HSD17B7 Antibody, Biotin conjugated

Product Code	CSB-PA010776LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56937
Immunogen	Recombinant Human 3-keto-steroid reductase protein (128-198AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	3-keto-steroid reductase (EC 1.1.1.270) (17-beta-hydroxysteroid dehydrogenase 7) (17-beta-HSD 7) (Estradiol 17-beta-dehydrogenase 7) (EC 1.1.1.62) (Short chain dehydrogenase/reductase family 37C member 1), HSD17B7, SDR37C1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HSD17B7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.